

Excellence Weighing Solutions



ErgoClips Unique Safe Dosing

ErgoClips allow easy weighing directly into tare containers



Don't waste your time
Say "Goodbye" to weighing paper. Its a thing of the past.



Perfect Positioning
ErgoClips guarantee secure placement of different tare containers.



Safe Dosing
Dosing is fast and safe thanks to the adjustable MinWeigh Door.




Accurate Weighing
With ErgoClips you can rely on complete sample transfer.

METTLER TOLEDO

ErgoClips Overview





Micro
Balances
XP26
and XP56




Part	ErgoClip Basket micro	ErgoClip Flask micro	MinWeigh Door micro	Stand for ErgoClips micro
Size	Ø max. 23 mm	1-10 ml	Inner draft shield with a small adjustable opening	Holder to store up to 3 ErgoClips
Tare containers (examples)	Vials, test tubes, centrifuge tubes	Volumetric flasks		
Order No.	11107889	11107879	11107869	11140175

XS/XP Analytical Balances



Part	ErgoClip Weighing Boat	ErgoClip Titration Basket	ErgoClip Basket small	ErgoClip Tube	ErgoClip Filter Holder	ErgoClip Quantos and vial adapters
Size	LxW: 60x45 mm	Ø max. 50 mm	LxW: 32x32 mm	Ø 8 / 12 mm	Ø 110-150 mm	Ø 10.3-21 mm
Tare containers (examples)	Weighing boats	Titration beakers, test tubes	Test tubes, centrifuge tubes or small plastic tubes	PCR tubes, Eppendorf® tubes, reaction tubes	Filters	HPLC tubes, vials
Order No.	11106748	11106883	11106747	11106784	11140185	11141570



Part	ErgoClip Round Bottom Flask	ErgoClip Flask	ErgoClip Small Flask	MinWeigh Door	Stand for ErgoClips
Size	Ø min. 46 mm	20-100 ml	1-10 ml	Draft shield with a small adjustable opening	Holder to store up to 3 ErgoClips
Tare containers (examples)	Round bottom flasks, beakers	Volumetric flasks	Volumetric flasks		
Order No.	11106746	11106764	11140180	11106749	11140170

Mettler-Toledo AG
Laboratory & Weighing Technologies
CH-8606 Greifensee, Switzerland
Phone +41 44 944 22 11
Fax +41 44 944 30 60

Subject to technical changes
© 03/2009 Mettler-Toledo AG
Printed in Switzerland 11796297
Global MarCom Greifensee

www.mt.com/ergoclips

For more information